# **Handling SAM files**

Estudios in silico en Biomedicina (Máster en Bioinformática, Universidad de Valencia)

David Montaner

2014-10-21

# **SAM Format**

Header section: @

#### Alignment section:

- QNAME: read name
- FLAG: bitwise FLAG
- 3 RNAME: reference sequence name (chromosome)
- POS: 1-based leftmost mapping position (where the read starts mapping)
- MAPQ: mapping quality
- CIGAR: CIGAR string
- RNEXT: name of the mate/next read (for paired ends)
- PNEXT: Position of the mate/next read (BP O )
- TLEN: observed Template LENgth
- SEQ: read sequence
- QUAL: read quality (phred-scale +33)

## **BAM** format

#### **SAM**

SAM: Sequence Alignment/Map format.

Is a TAB-delimited **text** format storing the alignment information.

#### **BAM**

A **binary** and *compressed* version of the file is easier to handle.

Generally **sorted** by chromosome position.

Supplementary **index** files may be created for certain purposes (quick access)

(quick access)

#### **SAMtools**

Program: samtools (Tools for alignments in the SAM format)

Version: 0.1.19-96b5f2294a

Usage: samtools <command> [options]

Command: view SAM<->BAM conversion

sort sort alignment file

mpileup multi-way pileup
depth compute the depth
faidx index/extract FASTA

tview text alignment viewer

index index alignment

idxstats BAM index stats (r595 or later)

fixmate fix mate information

flagstat simple stats

#### **SAMtools**

Program: samtools (Tools for alignments in the SAM format)

Version: 0.1.19-96b5f2294a

Usage: samtools <command> [options]

Command: calmd recalculate MD/NM tags and '=' bases

merge merge sorted alignments rmdup remove PCR duplicates

reheader replace BAM header

cat concatenate BAMs

bedcov read depth per BED region

targetcut cut fosmid regions (for fosmid pool or

phase phase heterozygotes

bamshuf shuffle and group alignments by name

## **SAMtools Commands**

```
Usage: samtools sort [options] <in.bam> <out.prefix>

Options: -n sort by read name
-f use <out.prefix> as full file name instead
-o final output to stdout
-l INT compression level, from 0 to 9 [-1]
-@ INT number of sorting and compression thread
-m INT max memory per thread; suffix K/M/G reco
```

# **SAMtools First Commands**

- view SAM- BAM conversion
- sort sort alignment file
- index index alignment

- reheader replace BAM header
- cat concatenate BAMs
- merge merge sorted alignments

- mpileup multi-way pileup
- depth compute the depth

## **SAMtools more Commands**

- faidx index/extract FASTA
- tview text alignment viewer
- idxstats BAM index stats (r595 or later)
- fixmate fix mate information
- flagstat simple stats
- calmd recalculate MD/NM tags and '=' bases
- rmdup remove PCR duplicates
- bedcov read depth per BED region
- targetcut cut fosmid regions (for fosmid pool only)
- phase phase heterozygotes
- bamshuf shuffle and group alignments by name